



Safe and Quick Lifting in a Horizontal Position.

Transferring patients in a horizontal position, for instance, in connection with surgery, radiology or intensive care, is often a challenge. It should be quick, simple and ergonomic for caregivers and comfortable and safe for the patient.

LikoStretch 1900.™ meets all of these requirements. The stretcher is connected to a patient lift and the patient can therefore be transferred in a horizontal position, e.g., between bed and gurney, without strenuous, manual lifting for the caregivers. Also, patient care, such as hygiene and sheet changing, is facilitated if the patient can be lifted with a stretcher rather than manually.

Furthermore, the open design allows for other equipment, such as tubes for respiratory aid, to be connected during the entire transfer operation.

Important Features

- Can be folded for smart storage
- Stretch leveller for optimal balancing
- Flat underlying lifting surface
- Simple strap lock system
- Full or 2-piece lift sheet as accessory

Fits Several Lifts

LikoStretch 1900™ can be used together with the overhead lifts Likorall™ and Multirall™, and the mobile lifts Viking™ L, Viking™ XL and Golvo™ ...

Accessories

Quick-Release Hook

If the lift and LikoStretch 1900™ are equipped with a Quick-Release Hook, the stretcher can be easily assembled when necessary, without the use of tools. LikoStretch 1900™ can also be combined with LikoScale™ to weigh patients in connection with lifting.

Stretch Leveller

With the included Liko™ Stretch Leveller™ you can easily adjust the lifting so that the patient lies either horizontally or slightly higher at the head end.



The lift sheet is attached to the stretcher by inserting the straps in the self-locking feature.



The fact that LikoStretch 1900 weighs only 19 kg (42 lbs) makes it easy to handle for one person. It can also be folded for convenient storage.

Lift Sheets

You can choose between a full or a 2-piece lift sheet. The 2-piece lift sheet can be especially suitable in connection with gynaecological operations, for instance, as the lower part can be taken away if necessary.

- Max. load:** 250 kg
- Length:** 190 cm (folded 70 cm)
- Width:** 90 cm (folded 90 cm)
- Height:** 60 cm (folded 70 cm)
- Weight:** 19 kg

Product Number:

- 3156051 LikoStretch 1900 incl. Stretch Leveller
- 3685705 LiftSheet 1900, full
- 3685704 LiftSheet 1900, 2-piece
- 3685905 Solo LiftSheet 1900, set of 10 pcs



LikoStretch 1900™ complies with the requirements for Medical Device Directive Class I products (MDD 93/42/EEC).

Design and Quality by Liko in Sweden.

Liko is quality certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485. Liko is also certified in accordance with environmental standard ISO 14001.

Hill-Rom and design and Hill-Rom word mark is a trademark or registered trademark of Hill-Rom Services Inc. Liko and design and Liko word mark is a trademark or registered trademark of Liko R&D AB. LikoStretch 1900 and Liko Stretch Leveller are trademarks or registered trademarks of Liko R&D AB.



Enhancing Outcomes for Patients and Their Caregivers.™

France _____ +33 (0)2 97 50 92 12
 United Kingdom _____ +44 (0)1530 411000
 Deutschland _____ +49 (0)211 16450-0
 Nederland _____ +31 (0)347 32 35 32
 Italia _____ + 39 (0)2 950541
 Suisse/Schweiz _____ +41 (0)21 706 21 30
 (deutschsprachig) _____ +41 (0)21 706 21 38

Österreich _____ +43 (0)2243 28550
 Ireland _____ +353 (0)1 413 6005
 Iberia _____ +34 (0)93 6856000
 Nordic Countries _____ +46 (0)20 78 10 30
 Export _____ Contact your local Hill-Rom distributor or contact your Area Manager via website or call _____ +1 812 934 8173

www.hill-rom.com
www.liko.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2010 Hill-Rom SARL, ALL RIGHTS RESERVED