



For a Wide Range of Use in Healthcare.

MultiStrap™ is a lifting accessory from Liko™ which completes the other lifting and repositioning aids in Liko's wide range. MultiStrap™ has been developed to meet a variety of needs in care and can be used to hold arms, legs and other parts of the body, for example, during treatment or washing.

Together with a patient lift, MultiStrap™ facilitates many heavy lifting tasks such as wound care (e.g. burn injuries), hygiene procedures, treatment (e.g. catheterizing or dressing wounds), examinations, placement of bedpans and dressing/undressing. The MultiStraps can also be used to turn a patient to the side.

MultiStrap™ can be used individually or in pairs. For turning in bed, we recommend using it in pairs with one of the wider models.

Important Features

- Facilitates lifting/holding e.g. arms and legs during treatment/washing
- Three washable MultiStraps in different sizes
- Also available are MultiStraps in disposable fabric intended for short-term use by an individual patient

Different Models and Areas of Use.

Liko™ MultiStrap™ Standard

Width 10 cm (4 inches)

Liko™ MultiStrap™ Standard is usually the easiest to apply and allows better access to a larger area of skin.

It is especially suitable for lifting and holding an arm or leg.

Liko™ MultiStrap™ Wide

Width 25 cm (10 inches)

Liko™ MultiStrap™ Wide has a greater support area and distributes pressure better thanks to its reinforcement.

It can be used to lift a patient's leg or behind, or to turn a patient in bed.

Liko™ MultiStrap™ Ultra

Width 40 cm (16 inches)

Developed especially for bariatric patients but can also be used to turn lighter patients in bed.

Has reinforcement to distribute pressure more evenly.



A FixStrap 135 cm (54 in) is included with all Liko™ MultiStraps and is used to fix one end of the MultiStrap™ to the bed when turning a patient.

Read the relevant instruction guides for further information on areas of use and practical application.

You can also contact your closest Liko™ representative for information and advice.

Product Information

Product Designation

Liko™ MultiStrap™ Standard, incl. FixStrap 3695115
 Liko™ MultiStrap™ Wide, incl. FixStrap 3695116
 Liko™ MultiStrap™ Ultra, incl. FixStrap 3695117

Product No.	Maximum Load	Fabric	Width
3695115	200 kg (440 lbs)	Polyester	10 cm (4 in)
3695116	200 kg (440 lbs)	Polyester	25 cm (10 in)
3695117	200 kg (440 lbs)	Polyester	40 cm (16 in)

Accessories

FixStrap 135 cm (54 in)
 Spatula 600 IC

Product No.

3696001
 3684002

(can facilitate placement of Liko™ MultiStrap™ Standard and Wide)

MultiStrap™ is a class 1 medical device product.

Design and Quality by Liko in Sweden.

Liko is quality-certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485.

Liko is certified in accordance with environmental standard ISO 14001.



Hill-Rom is a trademark or registered trademark of Hill-Rom Services Inc. Liko is a trademark or registered trademark of Liko R&D AB. MultiStrap is a trademark or registered trademark of Liko R&D AB.

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Manufacturer:
 Liko AB
 Nedre vägen 100
 S-975 92 Luleå,
 Sweden
 Tel. +46 (0)920 - 47 47 00
 Fax +46 (0)920 - 47 47 01
 info@liko.se

To find your local Liko representative,
 please visit www.liko.com or call
 +46 (0)920 47 47 00

www.hill-rom.com
www.liko.com

Hill-Rom® reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom® makes is the express written warranty extended on the sale or rental of its products.

© 2011 Hill-Rom® SARL, ALL RIGHTS RESERVED